

# SEMTRAC SET 4.2

## SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

### FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- three piece catheter for difficult transfers
- 5 free mandrels, to preshape the guiding catheter, in every box of 25 pieces
- reinforced shaft of loading catheter for optimal control
- handle with luer lock for a more secure fitting to syringes



ORIGINAL PICTURE

### GYNO INFO SEMTRAC SET 4.2

Guiding catheter (A) length 158 mm ( $\pm 2$  mm)  
tube OD 2.55 mm ( $\pm 0.05$  mm)


Loading catheter (C) length 232 mm ( $\pm 2$  mm)  
tube OD 1.70 mm ( $\pm 0.05$  mm)  
total volume 0.23 ml  
total volume minus hub 0.081 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

Order number #2000 Set 4.2 



GMP Belgium Semtrac Set 4.2 04/2011



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine