

SEMTRAC SET W-180

SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- 5 free mandrels, to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter

ORIGINAL PICTURE

GYNO INFO SEMTRAC SET W-180

Guiding catheter (A) length 141 mm (\pm 2 mm)
tube OD 2.55 mm (\pm 0.05 mm)

Loading catheter (C) length 214 mm (\pm 2 mm)
tube OD 1.70 mm (\pm 0.05 mm)
total volume 0.23 ml
total volume minus hub 0.074 ml


Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

Order number

- Semtrac Set W-180 #2000 W-180 ☒
- Semtrac Set W-180 Trial #2000 W-180 T ☒



GMP Belgium Semtrac W-180 12/2012

