

SEMTRAC SET W-240

SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- 5 free mandrels, to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- trial catheter set available for cervix examination prior to embryo transfer consisting of a guiding catheter and a soft closed-end trial catheter

ORIGINAL PICTURE

GYNO INFO SEMTRAC SET W-240

Guiding catheter (A) length 203 mm (\pm 2 mm)
tube OD 2.55 mm (\pm 0.05 mm)

Loading catheter (C) length 277 mm (\pm 2 mm)
tube OD 1.70 mm (\pm 0.05 mm)
total volume 0.26 ml
total volume minus hub 0.098 ml

Sterile ETO sterilized


Packaging 75/shipping box

Quality control MEA and LAL

Order number

- Semtrac Set W-240 #2000 W-240 

- Semtrac Set W-240 Trial #2000 W-240 T 



GMP Belgium Semtrac W-240 12/2012


Practise better Medicine

 Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be