

SEMTRAC VISION C

SOFT ECHOGENIC EMBRYO TRANSFER CATHETER

The Semtrac Vision C is a soft embryo transfer catheter, especially designed for ultrasound guided placement.


FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- one piece economical catheter
- markings to indicate depth
- metal tube inside for better control
- handle with luer lock for a more secure fitting to syringes
- echogenic catheter for optimal visualisation during ultrasound guided transfers



ORIGINAL PICTURE

GYNO INFO SEMTRAC VISION C

Length	252 mm (± 2 mm)
Tube OD	1.70 mm (± 0.05 mm)
Total volume	0.20 ml
Total volume minus hub	0.065 ml
Sterile	ETO sterilized
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#2000 SV C 



GMP Belgium Semtrac Vision C 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine