

SEMTRAC VISION SET 5

SOFT ECHOGENIC EMBRYO TRANSFER CATHETER

The Semtrac Vision Set is a soft embryo transfer catheter set for difficult embryo transfers, especially designed for ultrasound guided placement.

FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- three piece catheter for difficult transfers
- 5 free mandrels, to preshape the guiding catheter, in every box of 25 pieces
- reinforced shaft of loading catheter for optimal control
- handle with luer lock for a more secure fitting to syringes
- echogenic loading catheter for optimal visualisation during ultrasound guided transfers



ORIGINAL PICTURE

GYNO INFO SEMTRAC VISION SET 5

Guiding catheter (A) length 178 mm (± 2 mm)
tube OD 2.55 mm (± 0.05 mm)


Loading catheter (C) length 252 mm (± 2 mm)
tube OD 1.70 mm (± 0.05 mm)
total volume 0.22 ml
total volume minus hub 0.089 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

Order number #2000 SV Set 5



GMP Belgium Semtrac Vision Set 5 04/2011

