

# SEMTRAC VISION SET W-180

## SOFT ECHOGENIC EMBRYO TRANSFER CATHETER

The Semtrac Vision Set is a soft embryo transfer catheter set for difficult embryo transfers, especially designed for ultrasound guided placement.

### FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- 5 free mandrels, to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- echogenic loading catheter for optimal visualisation during ultrasound guided transfers



ORIGINAL PICTURE

### GYNO INFO SEMTRAC VISION SET W-180

Guiding catheter (A) length 141 mm ( $\pm$  2 mm)  
tube OD 2.55 mm ( $\pm$  0.05 mm)


Loading catheter (C) length 214 mm ( $\pm$  2 mm)  
tube OD 1.70 mm ( $\pm$  0.05 mm)  
total volume 0.23 ml  
total volume minus hub 0.074 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

Order number #2000 SV W-180 ☒



GMP Belgium Semtrac Vision W-180 04/2011

