

SEMTRAC VISION SET W-240

SOFT ECHOGENIC EMBRYO TRANSFER CATHETER

The Semtrac Vision Set is a soft embryo transfer catheter set for difficult embryo transfers, especially designed for ultrasound guided placement.

FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- 5 free mandrels, to preshape the guiding catheter, in every box of 25 pieces
- handle with luer lock for a more secure fitting to syringes
- echogenic loading catheter for optimal visualisation during ultrasound guided transfers



ORIGINAL PICTURE

GYNO INFO SEMTRAC VISION SET W-240

Guiding catheter (A) length 203 mm (\pm 2 mm)
tube OD 2.55 mm (\pm 0.05 mm)


Loading catheter (C) length 277 mm (\pm 2 mm)
tube OD 1.70 mm (\pm 0.05 mm)
total volume 0.26 ml
total volume minus hub 0.098 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

Order number #2000 Semtrac SV W-240 



GMP Belgium Semtrac Vision W-240 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine