

PREMIUM FAS

FOLLICLE ASPIRATION SET

The Premium FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.

FEATURES

- standard follicle aspiration needle with regular tubing
- sharp lancet shaped needle ensures a clean facet cutting
- etched ultrasound markings on tip which offer stable and excellent vision during ultrasound
- translucent tubing in the same size as needle for safe oocyte retrieval
- arrow marking on handle for bevel point indication
- handle shape facilitates rotation



ORIGINAL PICTURE

GYNO INFO PREMIUM FAS

Sterile sterilized by radiation

Packaging 50/shipping box
(standard configuration)

Quality control MEA and LAL

Order number **#4551 NS-AS1** ⓧ
#4551 NS-BS1 ⓧ
#4551 NS-CS1 ⓧ
#4551 NS-DS1 ⓧ



GMP Belgium Premium Fas Set 12/2012

GYNO INFO

PREMIUM SINGLE LUMEN

① Order number	Gauge	OD/ID (mm)	Length (mm)	Working length (mm)	Control
#4551 NS-AS1	17.5G	1.4/1.2	400	350	Flexible
#4551 NS-BS1	17.5G	1.4/1.2	340	290	Flexible
#4551 NS-CS1	17.5G	1.4/1.0	400	350	Rigid
#4551 NS-DS1	17.5G	1.4/1.0	340	290	Rigid

Order number	Gauge	OD/ID (inch)	Length (inch)	Working length (inch)	Control
#4551 NS-AS1	17.5G	0.055/0.047	15.75	13.779	Flexible
#4551 NS-BS1	17.5G	0.055/0.047	13.39	11.417	Flexible
#4551 NS-CS1	17.5G	0.055/0.039	15.75	13.779	Rigid
#4551 NS-DS1	17.5G	0.055/0.039	13.39	11.417	Rigid

- ② Tube from needle to cork : Regular tubing 100 cm / 39.369 inch
- ③ Size silicon cork : Small 14 x 18 x 20 mm / 0.551 x 0.709 x 0.787 inch
- ④ Tube from cork to pump : Regular tubing 100 cm / 39.369 inch

