

SUPERIOR FAS

CUSTOM-MADE FOLLICLE ASPIRATION SET

The Superior FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.


FEATURES



- sharp lancet or back cut bevel shaped needle ensures a clean facet cutting
- dual laser or embossed echogenic marked tip offers stable and excellent vision during ultrasound
- various needle sizes for ovaries with complicated and/or limited access
- translucent tubing in the same size as needle for safe oocyte retrieval
- arrow marking on handle for bevel point indication
- handle shape facilitates rotation
- custom-made



ORIGINAL PICTURE

GYNO INFO SUPERIOR FAS

Sterile	sterilized by radiation
Packaging	50/shipping box (standard configuration)
Quality control	MEA and LAL
Order number	#4551 



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine