

# GYNO INFO

## CUSTOM-MADE SUPERIOR FLUSH NEEDLES

### DOUBLE LUMEN FLUSH

Gauge	OD/ID (mm)	Length (mm)	Working length (mm)	Control
16G	outer lumen 1.65/1.40 inner lumen 1.20/1.00	352	300	Rigid
17G	outer lumen 1.50/1.20 inner lumen 1.00/0.80	352	300	Rigid

Gauge	OD/ID (inch)	Length (inch)	Working length (inch)	Control
16G	outer lumen 0.065/0.055 inner lumen 0.047/0.039	13.858	11.811	Rigid
17G	outer lumen 0.059/0.047 inner lumen 0.039/0.032	13.858	11.811	Rigid

### SINGLE LUMEN FLUSH

Gauge	OD/ID (mm)	Length (mm)	Working length (mm)	Control
17.5G	1.40/1.20	352	300	Flexible
17.5G	1.40/1.00	352	300	Rigid

Gauge	OD/ID (inch)	Length (inch)	Working length (inch)	Control
17.5G	0.055/0.047	13.858	11.811	Flexible
17.5G	0.055/0.039	13.858	11.811	Rigid

