

DETECTA

CERVICAL SCREENING DEVICE DESIGNED TO OBTAIN OPTIMAL SAMPLES

The Detecta is a simple but efficient cervical sampling device.

FEATURES

- allows in-office sampling
- systematic brushing of the transitional zone
- helps to detect neoplasia at an early stage
- easy to handle



ORIGINAL PICTURE

GYNO INFO DETECTA

Length	194 mm
Packaging	1250/shipping box
Order number	#4163 ☒



GMP Belgium Detecta 06/09



MANUAL DETECTA

- gently insert the Detecta into the cervix
- rotate only once to obtain the sample
- gently withdraw the Detecta without rotating
- smear sample on a slide with one rotation of the Detecta
- fixate the slide immediately
- sample is ready for examination

