

# STRIP-ID125-275

## MICROTOOL FOR DENUDATION OF EMBRYOS AND OOCYTES

The denuding pipette Strip-ID is a single-use laboratory microtool used for denuding of embryos and oocytes

### FEATURES











- superior quality that ensures smooth edges for precise denudation
- strict individual microscopic quality inspection
- flexible, unbreakable design that will not scratch dishes
- convenient multi-packaging that allows easy removal of a single pipette
- various sizes available to suit your needs
- compatible design to fit most commonly used pipette holders



ORIGINAL PICTURE

ORIG. PIC. 200X  
ENLARGED

### GYNO INFO STRIP-ID

Length	90 mm ( $\pm$ 3 mm)	
Inner diameter $\mu$ m	125 $\mu$ m (+ 3 $\mu$ m) 130 $\mu$ m (+ 3 $\mu$ m) 135 $\mu$ m (+ 3 $\mu$ m) 140 $\mu$ m (+ 3 $\mu$ m) 145 $\mu$ m (+ 3 $\mu$ m)	150 $\mu$ m (+ 3 $\mu$ m) 170 $\mu$ m (+ 3 $\mu$ m) 175 $\mu$ m (+ 3 $\mu$ m) 200 $\mu$ m (+ 3 $\mu$ m) 275 $\mu$ m (+ 3 $\mu$ m)
Sterile	sterilized by radiation	
Packaging	10 tips / tube 5 tubes / inner box 25 tubes of 10 tips / shipping box	
Quality control	MEA and LAL	
Order number	#3000 Strip-ID125/10  (125 $\mu$ m) #3000 Strip-ID130/10  (130 $\mu$ m) #3000 Strip-ID135/10  (135 $\mu$ m) #3000 Strip-ID140/10  (140 $\mu$ m) #3000 Strip-ID145/10  (145 $\mu$ m) #3000 Strip-ID150/10  (150 $\mu$ m) #3000 Strip-ID170/10  (170 $\mu$ m) #3000 Strip-ID175/10  (175 $\mu$ m) #3000 Strip-ID200/10  (200 $\mu$ m) #3000 Strip-ID275/10  (275 $\mu$ m)	

GMP Belgium Denuding Pipette Strip-ID 06/2013



Gynetics Medical Products n.v.  
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium  
tel. +32 11 64 58 72 / fax +32 11 64 06 77  
info@gynetics.be / www.gynetics.be

**GYNETICS**<sup>®</sup>  
Practise better Medicine