

EMTRAC DELPHIN

EMBRYO TRANSFER CATHETER WITH POSITIONER AND SOFT UTERINE CATHETER

The combination of a soft flexible intrauterine catheter and a solid cervix catheter creates an universal catheter for easy insertion in the most difficult cervical tract.


FEATURES


- positioner to prevent damage to fundus
- microscopic quality inspection of each tip
- markings on the outer cannula indicating uterine depth
- extra soft uterine catheter to ensure safe insertion through cervix
- handle with luer lock for a more secure fitting to syringes



ORIGINAL PICTURE

GYNO INFO EMTRAC DELPHIN

Length	210 mm (\pm 3 mm)
Tube OD	2.00 mm (\pm 0.20 mm)
Total volume	0.23 ml
Total volume minus hub	0.072 ml
Sterile	ETO sterilized
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4219 Delphin 



GMP Belgium Emtrac Delphin 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine